


Deep anesthesia causes poor perfusion, hypoxia and hypoventilation which can result in death. Anesthetic depth is difficult to assess based on physical signs. Heart and respiratory rate changes don't correlate with changes in Anesthetic depth, although monitoring respiratory rate is important because apnoea can lead to hypoxia. Better indicators are cardiac function (blood flow or pressure) and respiratory function (CO₂ clearance) which are depressed in a dose-dependent

Respiratory Function


Adequacy of ventilation depends on the respiratory rate AND depth (tidal volume) which determine minute ventilation and therefore arterial CO₂. Normal minute ventilation (about 200 ml/kg/min for dogs and cats) in conscious animals with normal lungs results in an arterial and therefore alveolar CO₂ partial pressure of 35 to 45 mm Hg.



Ventilatory function (minute ventilation) decreases and therefore CO₂ increases in a

dose-dependent fashion with increasing Anesthetic depth.

Comparison of Ventilatory Effects of Halothane, Isoflurane, and Methoxyflurane in Dogs Spontaneous Ventilation



End Tidal CO₂ sampling

At the end of expiration, gas sampled from the trachea should be alveolar gas and the CO₂ level will equate to the alveolar and therefore arterial CO₂ in dogs and cats that don't have lung disease.

End Tidal CO₂ (ETCO₂) monitors measure the CO₂ partial pressure using proportional absorption of infrared light. The monitors continuously sample gas at the endotracheal tube connector. The CO₂ concentration can be measured directly via a sensor in-line, called mainstream sampling, as shown below.




Alternatively the gas sample may be continuously drawn via a long sample line (sidestream sampling) up to the monitor where the CO₂ concentration is measured



ETCO₂ pitfalls and interpretation

A continuous CO₂ wave form can be generated which will demonstrate the three phases of expiration during each respiratory cycle.




- A = dead space gas, no CO₂
- B = mixed gas, rising CO₂
- C = alveolar gas, CO₂ plateau
- D = inspiration no CO₂

The monitor will continuously read the CO₂ level (in % or mm Hg) and display the peak CO₂ value (ETCO₂) and respiratory rate.


The plateau observed at the end of the CO₂ expiratory waveform (phase "C") ensures that the peak ETCO₂ value is determined from an alveolar gas sample.

Small tidal volumes, such as with cats, may not allow adequate sampling resulting in waveforms without an "alveolar gas plateau" such as below.



In this case the ETCO₂ may be underestimating the true alveolar and arterial CO₂.

Cardiac Oscillations



During expiration in animals with large chests, heart beats can force small puffs of gas up the trachea causing "bumps" in phase "C" of the ETCO₂ waveform.

The monitor will read each "bump" as a new breath, so give an erroneous respiratory rate and variable ETCO₂ values. The highest CO₂ level is in the arterial blood so when cardiac oscillations are present, observe the waveform and note the highest ETCO₂ reading which will be closest to the arterial CO₂.

Gasping respiration such as observed in animals at light planes of anesthesia, will cause similar problems due to variations in the waveform.

Lung disease or atelectasis such as seen during **thoracotomies** cause an increase in "dead space" ventilation, which dilutes the alveolar gas CO₂ concentration resulting in erroneously low ETCO₂ estimations. ETCO₂ values measured during a thoracotomy should be regarded as erroneous unless there is an arterial blood gas measurement to prove correlation. This dead space error also occurs in **anesthetized horses** where ETCO₂ values can underestimate arterial CO₂ by 5 to 15 mm Hg.

No ETCO₂ reading can be due to:


- failed intubation
- disconnection between the ET tube and Y piece/ETCO₂ sampling point
- airway occlusion
- lack of pulmonary blood flow (e.g. cardiac arrest, pulmonary air embolis)

Low ETCO₂ can be due to:

- light anesthesia and pain which increases spontaneous ventilation
- excessive mechanical ventilation
- small tidal volume causing a sampling error (eg cats)
- increasing dead space such as cause by alveolar collapse during a thoracotomy

High ETCO₂ can be due to:

- deep anesthesia
 - increased dead space (eg long ET tubes)
 - inadequate removal of CO₂ by exhausted soda-lime
 - low fresh gas flows in non-rebreathing circuits
 - rebreathing expired gas such as caused by defective 1-way valves
- High CO₂ causes high HR, BP and injected mucous membranes**



	ETCO ₂ Values, mm Hg		
	Anesthesia, Spontaneous Ventilation		
	Conscious	Light GA	Deep GA
Dog	35 - 45	45 - 55	55 +
Cat	30 - 38	35 - 45	50 +
Foal 2 mo.	35 - 45	45 - 55	60 +
Horse	35 - 45	55 - 65	75 +